

Gram Scales

Model SRV320A

When it's time to look to the sky, do not let the SRV320A pass you by. This scale allows veterinarian and zoological specialists to gather weights on birds up to 5,000 grams.

The large base, along with its 5" perch, permit the bird to be comfortably placed for proper weighing.

The perch is attached firmly to the base, however it can be removed allowing use of the scale for smaller animals or as a daily nutritional scale.



DISPLAY RESOLUTION

5000 grams x 1 gram resolution.

DISPLAY UNITS

Grams or ounces.

5" PERCH

Allows larger birds room to perch.

HOLD

Freezes the weight on the display.

BATTERY OPERATED

4 "AA" batteries. Auto Power Down (APD) approximately 120 seconds.

OPTIONAL EQUIPMENT

AC power supply

Gram Scales

Model SRV320A

SR Scales[®]
by SR Instruments, Inc.

Quality. Reliability. Accuracy.

SYSTEM SPECIFICATIONS

MAX WEIGHT 5000g / 175 oz

DISPLAY RESOLUTION 1g / 0.1oz

DISPLAY UNITS Grams / Ounces

POWER 4 "AA" alkaline batteries

AUTO POWER DOWN (APD) Approximately 120 seconds to conserve battery life. A low battery indication will appear when it's time to replace the batteries.

SCALE WEIGHT 2 Lb (1 Kg)


WARRANTY 2 years



ACCESSORIES AND OPTIONS

DESCRIPTION	PART NUMBER
AC ADAPTER	FRADATH004

DIMENSIONS in [cm]



PACKAGING

