

Heavy Duty Scale

Model SRV942

When you need a scale that is durable and heavy duty, consider the SRV942. While used primarily in stationary veterinary practices, this scale also has a wide range of applications in the zoological field.

The durable steel frame and surface are coated with SR's non-slip polyurea treatment and offer a 550 pound weighing capacity. *Stainless steel platform top available as option.

The rechargeable battery allows remote weighing when required. The SRV942 is equipped with a wall mount AC adaptor for stationary applications.

The scale is equipped with wheels and a fixed handle to allow for easy transport when required.

NON-SLIP PLATFORM

Durable and easy to clean polyurea coating over stainless with aggregate mixed in to create a surface animals are comfortable walking on. Optional Stainless top available.

PLATFORM SIZE

22" x 36" platform

AUTO-POWER SHUTDOWN

Automatically shuts down after two minutes of non use when not wall powered.

MULTIPLE POWER SOURCE OPTIONS

Rechargeable battery with AC/DC Adaptor



Heavy Duty Scale

Model SRV942



SR Scales[®]
by SR Instruments, Inc.

Quality. Reliability. Accuracy.


SYSTEM SPECIFICATIONS

MAX WEIGHT	550 Lb / 249 Kg
DISPLAY RESOLUTION	.2 lb / .1 kg
BUTTON FUNCTIONS	On/Zero/Off, Tare, Unit, and Hold
DISPLAY UNITS	Pounds / Kilograms
HOLD	Freezes the weight on the display
TARE	Allows objects to be placed on scale and zeroed out.
POWER	Rechargeable battery with DC Adaptor
SCALE SURFACE	SRV942 Polyurea surface w/ aggregate SRV942SS Stainless Steel
AUTO-POWER DOWN	Approximately two (2) minutes of non-use
PORTABILITY	Built-in handle and wheels
SCALE WEIGHT	Scale weighs 38 Lbs.
WARRANTY	Two (2) Years

DIMENSIONS in [cm]



PACKAGING



ACCESSORIES AND OPTIONS

DESCRIPTION	PART NUMBER
STAINLESS STEEL TOP	SRV942SS
REPLACEMENT CHARGER	SRV942CH

Specifications subject to change without notice.